


The only Medical Device of its kind  
for cervical cytology

Hygienic  
and  
Single-use  
bottle



# Colposol<sup>®</sup>


Solution 5%

Colposol<sup>®</sup> is a 5% w/v solution of glacial acetic acid which is used as a diagnostic adjunct in cervical cytology.

Legal Category: Class 1 Medical Device

Another trusted  
product from






**Colposol<sup>®</sup>** is the only Class 1 Medical Device of this formulation available for use as a diagnostic adjunct in cervical cytology.

**What is the purpose of Colposol<sup>®</sup>?**

**Colposol<sup>®</sup>** is used in the visual inspection of the cervix (the neck of the womb) when screening for cervical cancer. During the procedure the cervix is washed with **Colposol<sup>®</sup>** which removes any mucus and helps to identify any abnormal cells by staining them white, thereby differentiating them from the normal pink-coloured cells.

**Use Colposol<sup>®</sup> with confidence:** Single use bottles to minimise cross-contamination.

**Infection control is key for patient safety and quality of care.**



## Benefits of Colposol®:

- ✓ No need to decant from a larger vessel
- ✓ Practical single-use bottle to minimise cross-contamination and support infection control
- ✓ UK manufactured
- ✓ Stocked in all major wholesalers nationwide.
- ✓ Tested chemically and microbiologically
- ✓ Brought to you by Derma UK, trusted for over 20 years

Single use,  
for professional  
use in the  
hospital setting  
only



# Colposol®

Solution 5%

Colposol® is a 5% w/v solution of glacial acetic acid which is used as a diagnostic adjunct in cervical cytology.

Legal Category: Class 1 Medical Device

Another trusted  
product from



**Legal Category:**  
Class 1 Medical Device

**Pack Sizes:**  
30ml Bottle

**NHS Pricing:**  
£11.24

**For hospital use only:**  
For queries relating to  
**Colposol®**, please  
contact us on  
+44 (0)191 375 9020.



To find out more information about

**Colposol®**

Solution 5%

please visit [www.dermauk.co.uk](http://www.dermauk.co.uk)

Another trusted  
product from



**Adverse Effects**

Adverse events should be reported. Information about adverse event reporting can be found at [www.mhra.gov.uk/yellowcard](http://www.mhra.gov.uk/yellowcard). Adverse events should also be reported to Derma UK Ltd, UK on 0191 375 9020.